

CAPRICORN LABORATORY EQUIPMENT LIMITED
LEC MEDICAL REFRIGERATION RANGE – PRICES EFFECTIVE 1 SEPTEMBER 2009

MODEL	TEMPERATURE RANGE	CAPACITY (LITRES)	DIMENSIONS (MM) WxDxH	PRICE
-------	-------------------	-------------------	-----------------------	-------

PHARMACY REFRIGERATORS SOLID DOOR

PE102 (lockable)	+2°C to +8°C	41	379 x 500 x 702	£470
PE207 (lockable)	+2°C to +8°C	82	502 x 540 x 660	£610
PE507 (lockable)	+2°C to +8°C	153	558 x 600 x 863	£715
PE507DP (4 drawers)	+2°C to +8°C	153	558 x 600 x 863	£1250
PE907 (lockable)	+2°C to +8°C	286	595 x 600 x 1550	£940
PE907DP (8 drawers)	+2°C to +8°C	286	595 x 600 x 1550	£1855
PE357 (lockable)	+2°C to +8°C	354	595 x 600 x 1850	£1385
PE1607 (lockable)	+2°C to +8°C	444	685 x 910 x 1688	£2445

PHARMACY REFRIGERATORS GLASS DOOR

PG102 (lockable)	+2°C to +8°C	41	379 x 500 x 702	£560
PG207 (lockable)	+2°C to +8°C	82	502 x 508 x 660	£815
PG507 (lockable)	+2°C to +8°C	153	558 x 600 x 863	£925
PG907 (lockable)	+2°C to +8°C	286	595 x 600 x 1550	£1225
PG1607 (lockable)	+2°C to +8°C	444	685 x 910 x 1688	£2660

LABORATORY REFRIGERATORS

LR207 (lockable)	+2°C to +8°C Spark-free	82	502 x 540 x 660	£375
LR307 (lockable)	+2°C to +8°C	107	495 x 580 x 830	£435
LR507 (lockable)	+2°C to +8°C	153	558 x 600 x 863	£280
LR907 (lockable)	+2°C to +8°C	286	595 x 600 x 1550	£385
LR1607 (lockable)	+2°C to +8°C Spark-free	475	685 x 910 x 1688	£2340

LABORATORY FREEZERS

ISU27 (lockable)	to -20° Spark-free	50	502 x 540 x 660	£480
ISU37 (lockable)	-18°C to -25°C	102	495 x 580 x 830	£460
ISU57 (lockable)	-18°C to -25°C	116	558 x 600 x 863	£335
ISU97 (lockable)	-18°C to -25°C	230	595 x 600 x 1550	£505

LABORATORY FRIDGE/FREEZERS

ISR27	+2°C to +8°C/minimum to -9°C Spark-free	78/6	502 x 540 x 660	£375
ISR57 (lockable)	+2°C to +8°C/minimum to -18°C	112/18	558 x 600 x 863	£280
IST47 (lockable)	+2°C to +8°C/minimum to -25°C	157/89	560 x 600 x 1572	£555

LOW & ULTRA LOW TEMPERATURE REFRIGERATION

CLT111 (lockable)	-20°C to -45°C	113	730 x 850 x 1080	POA
CLT112 (lockable)	-40°C to -85°C	113	730 x 850 x 1080	POA
ULT341 (lockable)	-20°C to -45°C	322	765 x 900 x 1985	POA
ULT342 (lockable)	-40°C to -85°C	322	765 x 900 x 1985	POA
ULT531 (lockable)	-20°C to -45°C	500	1015 x 900 x 1985	POA
ULT532 (lockable)	-40°C to -85°C	500	1015 x 900 x 1985	POA

WARD REFRIGERATION

WR207 (lockable)	+1°C to +5°C	82	502 x 540 x 660	£425
WR507 (lockable)	+1°C to +5°C	153	558 x 600 x 863	£425
WR907 (lockable)	+1°C to +5°C	286	595 x 600 x 1550	£525

EXTRAS – TEMPERATURE MONITORING

USB Data Logger	Measurement Range -35°C to +80°C; up to 16382 Readings; Long Life Battery			£55
-----------------	---	--	--	-----

For further information, brochures or pdf files, please contact us by telephone on **01425 650088** or e-mail info@capricorn-labs.co.uk Alternatively, please download the information from our website at : www.capricorn-labs.co.uk

All prices are subject to VAT at the current rate. Delivery to ground floor Goods Inwards is £25 + VAT. Delivery to other locations, to include unpacking and removal of packaging from the premises, is £50 + VAT. Delivery prices for the Low & Ultra Low Temperature Refrigeration Range will be quoted on request.

In line with Glen Dimplex's policy of continual improvement, specifications are subject to change without notice. Capricorn reserves the right to alter prices without prior notice.